



THE LEADERS IN RECREATION, FITNESS AND ABILITY!

9889 Garrymore Ln
Missoula, MT 59808
888-687-3552 | +1-406-549-0769
www.aquacreek.com



Lift Cover

FOR:

PRO POOL-AT

RANGER-AT

AMBASSADOR

ADMIRAL

PATHFINDER

PORTABLE PRO POOL-AT

PATRIOT-AT

PART #: F-120PPC (Blue)

PART #: F-120PPC-T (Tan)

PART #: F-120PPC-G (Gray)

MANDATORY

LEAVE THIS MANUAL WITH LIFT OWNER

Read and follow all instructions.

Lift safety can only be ensured if the lift is installed and operated according to these instructions.

- NEVER permit children to play on or around the lift
- Do not allow children to use the lift without adult supervision
- NEVER apply direct water pressure to the electronics
- NEVER use the lift with a dry pool

READ CAREFULLY

Check entire box and all packing materials for parts. Before beginning assembly, read the instructions and identify parts using the figures and parts listed in this document.

It is critical that all parts be carefully inspected prior to installation. If any damage occurred in transit, Aqua Creek Products, LLC must be notified within three days of receipt of unit.

Proper installation cannot be overstressed, as an improper installation voids Aqua Creek's warranty and may affect the safety of the user.

Fixed Lift Cover Installation

For the Pro Pool-AT, Ranger-AT, Ambassador, Admiral and Pathfinder lifts, the cover is slid over the lift, then pull the draw string tight. The cover has grommet holes to allow for locking. PHOTO 'A' shows the cover installed over a fixed lift.



COVER IS OVER THE LIFT; USE THE
DRAW-STRING TO TIGHTEN

PHOTO 'A'

Portable Lift Cover Installation

For the Patriot-AT and Portable Pro Pool-AT lifts, the unit needs to be in the TRANSPORT mode, with the lift on the cart, and the cover fits over the pool-lift portion of the unit. PHOTO 'B' shows the cover installed over a portable lift.



COVER IS OVER THE LIFT AND
HANDLE; USE THE DRAW-STRING
TO TIGHTEN

PHOTO 'B'